Ref. : IFU_H01_05.23

NOTICE D’INSTRUCTIONS

Tiges CERAFIT HAC

Français

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CERAFIT HAC stems

English

Instructions for the following implants:

- CERAFIT MULTICONES HAC stems (ref. 4307 to 4318)
- CERAFIT R-MIS HAC stems (ref. 12007 to 12018)
- CERAFIT R-MIS HAC lateralized stems (ref. 12207 to 12216)

Graphic symbols used for the labelling of these devices:

Symbol	Title of the symbol	Symbol	Title of the symbol
	Manufacturer		Do not reuse
	Use-by date		Terminally sterilised device by irradiation.
	Batch code		Please note that not all the information needed to use the device is on the label, so it is necessary to consult the instruction manual.
	Catalogue reference of the device		

Before using a CERAVER® hip joint prosthesis, the surgeon must take note of the following recommendations:

1 - GENERAL

- The manufacturer and supplier of a CERAVER® hip prosthesis guarantee the quality of manufacture and materials. A prosthesis that has already been used must not be reimplanted. The assembly instructions for the individual components as given in the sales documentation must be observed. It is strictly forbidden to use CERAVER® hip prosthesis components together with components of another origin.
- Several criteria are of decisive importance for the success of an operation with a CERAVER® prosthesis:
 - Appropriate selection of the patient for the procedure.
 - A perfect surgical technique and correct use of the surgical instruments specially developed by CERAVER®.
 - Strict aseptic conditions.
- The surgeon is responsible for complications resulting from incorrect indications and faulty surgical technique or aseptis, which cannot be attributed to the manufacturer or supplier of CERAVER® products.
- The surgeon must be familiar with all the problems involved in the implantation of a joint replacement. The surgical technique for the implantation of a CERAVER® hip prosthesis can be taught in a surgical centre where the surgeons have acquired a perfect command of it.
- Patients receiving a hip joint replacement should be advised that their weight and activity level may affect the longevity of the prosthesis.

2 - DESCRIPTION OF THE IMPLANTS

Cementless femoral stems are manufactured using the latest technology. Their mechanical and functional characteristics can only be guaranteed if the operating technique and the instructions in this manual are followed.

These femoral components are designed for cementless implantation.

- The femoral stems are made of TiAl6V4/ELI titanium alloy (ISO 5832-3, ASTM F136).
- Femoral components for cementless fixation have a rough surface coated with hydroxyapatite (HA) (standards: ISO 13779-2, ISO 13779-3 and ISO 13779-4).

3 - INDICATIONS

- Hip replacement may be considered when all other surgical options have been carefully weighed and considered less appropriate.
- A joint prosthesis, even if successfully installed, always has inferior qualities to those of a natural joint. A joint prosthesis is an advantageous replacement for a pathological joint because it eliminates pain, restores good joint mobility and ensures that the operated limb is well supported.
- Every joint prosthesis is always subject to a degree of wear, which cannot be avoided. A prosthesis that was initially stable can become loosened over time. Wear and loosening of the prosthesis may lead the surgeon to re-operate.
- Infection of a prosthesis usually leads to its removal and sometimes to disabling sequelae.

Planned applications:

Hip arthroplasty is the replacement of a hip joint affected by various conditions that cannot be treated by other therapies or methods. The use of cementless femoral stems is aimed at restoring joint mobility, joint stability (absence of dislocation) and pain reduction (predicted performances). These performances are measured through a functional score, which is improved after surgery.

Hip joint replacement may be indicated in the following cases:

- Primary coxarthrosis.
- Secondary coxarthrosis (sequelae of trauma, infection, or necrosis) or favoured by congenital deformities (dysplasia) of the joint,
- Osteonecrosis of the femoral head,
- Non-infectious arthritis,
- Fracture of the proximal femur.

Based on the above considerations, a joint replacement is recommended when the following conditions are met:

- The anatomical shapes and quality of the bone structures in the patients must allow for the anchoring of a prosthesis.
- Patients should be informed by the surgeon of the risks of implanting a prosthesis and give their consent to the proposed procedure.

The following circumstances require special precautions:

- Overweight
- Excessive loading to be expected due to occupation or sports practices.
- Epilepsy or other factors leading to repeated accidents with an increased risk of fractures.
- Significant osteoporosis and/or osteomalacia.
- Old infectious diseases with potential for progression, with possible joint manifestations.
- Significant deformity of the joint to be replaced, suggesting difficulties in anchoring the prosthesis.
- Reduced bone support capacity as a result of tumour.
- Intellectual inability of patients to understand and follow the surgeon's instructions.
- Alcohol or drug abuse;
- Metabolic disorders or general drug treatments, leading to progressive deterioration of the bone support of the prosthesis.
- Conditions in other joints (knee, ankle,...).
- Overall weakened immune status of the patient.
- Unbalanced diet;
- Patients with small bones.

4 - CONTRAINDICATIONS

The following circumstances generally constitute an absolute contraindication to the placement of a joint prosthesis:

- Allergy to any component of the prostheses.
- Active infection of the joint or periarticular area in question, including severe neuroarthropathy;
- Loss of bone substance upstream or downstream of the joint, which makes proper anchoring of a cementless prosthesis impossible.
- Muscle loss, neuromuscular injury or vascular insufficiency of the affected limb, making the procedure unjustifiable;
- Charcot's disease;
- Patients under 50 years of age with joint damage that is likely to be successfully treated with conservative therapy;
- Pregnancy.

5 - IMPORTANT WARNINGS

5.1 Preoperative warnings

- The operation must be carefully planned. Depending on the results of the X-ray examination (tracings) and the preoperative planning, the surgeon will choose the appropriate implant.
- CERAVER® trays and instruments specially designed for the preparation of joint components should be used. The description of the surgical technique can never be exhaustive, nor can it include all the risks and complications that must be taken into account. Monographs on the individual procedures, technical instructions on our products and publications with more detailed information as well as video cassettes can be requested from CERAVER® or its distributor.
- The male taper of the femoral stem must match the female taper of the femoral head.
- The dimensions of the cone are indicated on the box and on the implant.

5.2 Intra-operative warning:

- Joint replacements should be handled with care. Scratches or irregularities caused by impact lead to excessive wear and tear and can cause complications.
- For the selection of the prosthesis size: only use CERAVER® trays, ancillary instruments and trial implants.
- Exposed prosthetic joint surfaces must not come into contact with any metal objects, such as instruments.
- It should be noted in particular that ceramic heads should be implanted on perfectly clean and dry cones.
- The protective cap of the cone should be removed prior to implantation of the femoral head.
- In terms of joint mobility, the worst case is the combination of a ø28 stainless steel femoral head +8 or +12 neck and an antiluxation core set in a Cerafit cup (regardless of the femoral stem). The lip core is positioned to artificially create a superior-lateral heel -generally stop and to limit the risk of dislocation. In this case, only the flexion/extension mobility is reduced by 18% compared to the minimum normative expectations in force (82.3° instead of the 100° required by NF EN ISO 21535). Similarly, when combining a ø32 stainless steel head with a +8 or +12 neck and an antiluxation core in the same positions, only the mobility in flexion/extension is reduced by 2% compared to the minimum normative expectations in force (98.1° instead of 100° required by NF EN ISO 21535). Mobility is within normative expectations when the core lip is positioned to create a posterior stop. This choice of implant combination and its impact should be taken into account before implantation.
- In terms of the mechanical fatigue strength of the stem neck, the most unfavourable case is the combination of a ø28 or ø32 stainless steel femoral head, most +12 (references 2262 and 2264) with a CERAFIT R-MIS lateralised stem size 16 (references 12416 and 12216). This configuration does not meet the requirements of the current NF ISO 7206-6 standard (the phenomenon observed is fatigue failure of the neck in less than 10 million cycles). This configuration is therefore prohibited, as it exposes the patient to a proven risk of fatigue failure of the neck of the femoral stem.

Warnings related to the anchoring technique

The anchoring is crucial to ensure reliable fixation of the implant. The following defects can lead to loosening and complications:

- Excessive thinning of the bone cortex during the preparation of the bone bed.
- inappropriate selection of the size of the prosthesis.
- Excessive force during the placement of the prosthesis fixation, resulting in burst fractures or bone splitting.

6 - POSTOPERATIVE PRECAUTIONS

- Every patient with a CERAVER® joint prosthesis must undergo a post-operative check-up by an orthopaedic surgeon. This is necessary in order to detect signs of wear or loosening of the prosthesis before functional manifestations occur. Early recognition of a complication allows effective action to be taken. In particular, an early reoperation has a better chance of success than a late one, especially when an advanced loosening has already led to significant resorption damage in the bone.
- The patient should report any changes in the operated limb to the surgeon.
- If an annual check-up is not possible, a control X-ray should be sent to the surgeon at least once a year for evaluation, as the X-ray study can detect early lesions.

6.4 Cerafit HAC cementless stems are non-ferromagnetic, have a geometry that is not likely to generate induced currents, and are fixed to the bone or implant, which reduces the risk of mobilisation. Their safety is particularly in terms of heating and migration of the implant, was assessed from literature data by comparison with data available on devices of similar

composition, shape and use. This evaluation concluded that the device is conditionally compatible for use in MRI scanners up to 3 Tesla. However, it is recommended that MRI examinations be avoided within 48 hours of implant placement and that the examiner be informed of recent implant placement if such an examination is required.

It should be noted that all devices with a high contrast with the surrounding biological environment may generate artefacts that should be considered for the perfect performance and interpretation of imaging examinations. To this end, it should be recommended that the patient with an implant informs the health professionals concerned (radiologists and radiology technicians) of the presence of an implant before the examination.

7 - COMPLICATIONS

The following complications have sometimes been observed and therefore require special attention:

Implant-related side effects:	Side effects of the procedure:
<ul style="list-style-type: none">• Fractures and cracks; • implant fracture; • pain; contractures • lower limb inequality of more than 1 cm; • sinking or migration of the stem; • Stem malposition; • impingement (pinching of the joint) • osteitic ossifications; • osteolysis, granuloma; foreign body reactions • dislocations; instability; • vascular lesions, haemorrhage • prosthesis infection • generation of noise in the prosthetic joint.	<ul style="list-style-type: none">• lung disorders, including pneumonia and atelectasis; • digestive complications paralytic ileus, inflammation and then necrosis of the intestine, which can lead to death; • delayed healing; • surgical wound hematoma; • mental confusion; • nerve palsy; superficial nerve damage, sciatica palsy (neurological injury); • torn tendon insertion (musculotendinous injury); • infection; • cardiovascular complications: atrial

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