Code article : 23000102  
Code Catalogue : H24  
Indice : 04.15

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## Français

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Cemented dual mobility cup

## English

**Instructions for use for the following implants:**  
- Cemented dual mobility cup (ref.: 14046, 14048, 14050, 14052, 14054, 14056, 14058, 14060, 14062)  
- PE insert for dual mobility cups (ref.: 10029, 10030, 10031, 10032, 10033, 10034, 10035, 10036, 10037)

Graphic symbols used for labelling these devices:

Symbol	Symbol title	Symbol	Symbol title
	Manufacturer		Do not re-use
	Expiration date		Device terminally sterilized by irradiation
	Batch code		Indicates a medical device that has been submitted to a sterilisation process
	Catalogue reference of the device		N.B.: Not all the information required for use of the device is indicated on the label; it is necessary to consult the instructions for use

**Prior to using a CERAVER<sup>®</sup> hip joint prosthesis, surgeons should familiarise themselves with the following recommendations:**

### 1- GENERAL POINTS

- 1.1** The manufacturer and the supplier of the CERAVER<sup>®</sup> hip prosthesis provide a warranty for the design and material quality of the product. A prosthesis which has already been used must not be implanted again. Instructions for assembling the different components must be respected, as defined in the commercial documentation. It is strictly forbidden to use CERAVER<sup>®</sup> hip replacement components together with implants of another origin.
- 1.2** The following factors are of prime importance for the success of an operation using a CERAVER<sup>®</sup> hip prosthesis:
  - Careful selection of a suitable patient for the operation in question.
  - Exact surgical technique and proper use of special CERAVER<sup>®</sup> instruments.
  - Extremely strict asepsis.
- 1.3** The surgeon is responsible for any complications that may result from an erroneous indication, incorrect surgical technique or inadequate aseptic precautions. Neither the manufacturer nor the supplier of the CERAVER<sup>®</sup> implants can be held responsible for these complications.
- 1.4** The surgeon must be completely familiar with all the aspects and problems of joint replacement surgery. The surgical procedures for the implantation of a CERAVER<sup>®</sup> prosthetic hip can be best learned at a clinic with many years' expertise in the use of these implants.
- 1.5** Patients receiving a hip joint prosthesis should be informed that their weight and level of activity may affect the longevity of the prosthesis.

### 2- DESCRIPTION OF IMPLANTS

CERAVER<sup>®</sup> acetabular components are manufactured using the latest technology. Their mechanical and functional characteristics can only be guaranteed if correct surgical technique and the present instructions are strictly observed. The cemented dual mobility cup is designed for implantation using orthopaedic cement.

- The inserts are made of ultra-high molecular weight polyethylene (UHMW PE) (NF ISO 5834-2 standard). The cemented dual mobility cup is made of wrought stainless steel with a high nitrogen content (NF ISO 5832-9 standard, ASTM F 1586 standard).

- The angular movements afforded by CERAVER<sup>®</sup> total hip replacements comply with the NF EN ISO 21535 standard.

- The angular movements afforded by CERAVER<sup>®</sup> total hip replacements are in compliance with NF EN ISO 21535.

### 3- INDICATIONS

- A joint prosthesis should be considered only after all other surgical methods of treatment have been carefully weighed up against each other and none have been judged to be more appropriate.
- Even if a patient has successfully implanted artificial joint is inferior to a sound natural joint. On the other hand, an artificial joint can be a highly beneficial substitute for a deformed or diseased joint, because it eliminates pain and is conducive to the restoration of good mobility and weight-bearing capacity.
- Every artificial joint is always subject to wear, which still remains a major problem that cannot be avoided. An initially stable prosthesis may become loose in the course of time. Wear and loosening of a prosthesis may make revision surgery necessary.
- Infection of a prosthesis will usually necessitate its removal or replacement in certain conditions and any such infection may leave debilitating residual effects.

### Planned applications:

Hip arthroplasty consists of fully or partially replacing a hip joint suffering from various affections that cannot be treated through other therapies or methods. It aims to restore all or part of the patient's mobility and stability while minimising pain or even eliminating it. The replacement of the hip joint may be indicated in the following cases :

- Deterioration of the joint resulting from primary or secondary osteoarthritis, traumatic arthritis, rheumatoid arthritis or dysplasia of the hip;
- Traumatic fracture of the head or neck of the femur;
- Avascular necrosis of the femoral head;
- Some cases of ankylosis;
- Etat après des interventions chirurgicales antérieures telles que la reconstruction d'une articulation (ostéotomie), une arthrodèse ou une hémiarthroplastie.

Given the above, it follows that a joint prosthesis is indicated in cases where the following conditions are fulfilled :

- The selected patient's joint is anatomically and structurally suited for the reception and implantation of a prosthesis;
- The patient must have been informed by the surgeon of the risks associated with the implantation of a joint prosthesis and has declared their consent to the operation.

Specific indications for the cemented dual mobility cup are :

- Risk of dislocation due to muscle deficiency, dysplasia, injury, or the physical or mental incapacity of the patient;
- Hip replacement revision;
- Rehabilitation following acetabular destruction enabling cementing of the cup in the acetabulum or in a metallic cage (e.g. a Kerboull-type acetabular reinforcement device, etc.).

### The following circumstances call for particular caution:

- Overweight patients;
- Excessive loading to be expected as a result of work or sport;
- Epilepsy or other factors favouring repeated accidents with increased risk of fracture;
- Severe osteoporosis and/or osteomalacia;
- Past history of still potentially active infectious diseases with potential joint manifestations;
- Severe deformation of the affected joint, which may render fixation of the implant more difficult;
- Fracture instability of the head or neck of the femur;
- Drug or alcohol addiction;
- Metabolic disorders or general medical treatments, leading to a progressive deterioration of bone support of the prosthesis;
- Other joint affections (knee, ankle, etc.);
- Patient's overall immune status weakened;
- Unbalanced diet.

### 4- CONTRAINDICATIONS

The following conditions are generally accepted as contraindications for the implantation of a joint prosthesis:

- Allergy to the implant materials of which the prosthetic components are made;
- Patients with a native cup with a diameter of less than 46 mm or greater than 62 mm;
- Active infection of the joint or surrounding area;
- Severe neuroarthropathy;
- Loss of bone substance above or below the joint, making correct anchoring of the cup or the metallic cage impossible;
- Significant factors may make joint prosthesis stability impossible, neuromuscular damage or vascular insufficiency of the limb affected, making the surgery unjustifiable;
- Charcot joint disease;
- Patients aged less than 50 years old and suffering from joint damage that could probably be treated successfully with a conservative treatment;
- Pregnancy;
- Tumours.

### 5- IMPORTANT WARNINGS

**5.1 Preoperative warnings**  
- The operation should be planned carefully. On the basis of the radiological templates and pre-operative planning, the surgeon will choose the appropriate implant.  
- The special CERAVER<sup>®</sup> templates and ancillary instrumentation should be used to prepare the fitting and setting of the articular components. A description of the operative technique can never be exhaustive and does not take into account all the risks and complications that should be considered.  
- CERAVER<sup>®</sup> PE dual mobility inserts must be used with a dual mobility cup of the same diameter.

#### 5.2 Intraoperative warnings

- Joint prostheses should be handled with care. Even slight scratches or impact dents on the articular surfaces may contribute to excessive wear and initiate undesirable complications.  
- To choose the right-sized implant: only CERAVER<sup>®</sup> templates, ancillary instruments and dummy models may be used.  
- Exposed articular surfaces must not come into contact with any metallic object (e.g. instruments).  
- Assembly of the insert and corresponding head must be performed in accordance with surgical technique using CERAVER<sup>®</sup> instruments.

#### 5.3 Warnings linked to the fixing technique

Proper fixation is of paramount importance for the firm and long-term seating of the implant. The following factors may make joint prosthesis stability impossible, neuromuscular damage or vascular insufficiency of the limb affected, making the surgery unjustifiable:

- Excessive weakening of the bone structure during the surgical preparation of the osseous bed;
- Infection of unsuitable prosthesis size;
- Inadequate pressure exerted on the implant, or mobilisation of the implant during the cement
- Inadequate application of force during insertion or positioning of the prosthesis, which may induce forced fractures or tears in the bone.

The cementing technique must also be prepared meticulously. It is particularly important to avoid the following errors or deficiencies:

- Non-compliance with the cement preparation instructions and, in particular, incorrect dosage or an excessively high initial temperature of the cement ingredients.
- Mixing in incorrect proportions or incorporating air bubbles when mixing the polymer and monomer.
- Careless cement application, with incorporation of blood and air.
- Inadequate pressure exerted on the implant, or mobilisation of the implant during the cement
- Inadequate application of force during insertion or positioning of the prosthesis, which may induce forced fractures or tears in the bone.

### 6- POSTOPERATIVE PRECAUTIONS

**6.1** All patients with an articulating CERAVER<sup>®</sup> prosthesis must be postoperatively monitored by an orthopaedic surgeon. Follow-up examinations are essential to detect signs of wear or prosthesis loosening before any resulting detrimental effects. Early detection of an impending complication allows the surgeon to initiate timely and effective countermeasures. Above all, a re-operation performed at the right time has a much better chance of success than a later surgical intervention, when advanced loosening has already resulted in serious resorption of the bone bed.

**6.2** Patients must be instructed to inform their surgeon of the slightest changes in the joint operated on.

**6.3** If annual clinical check-ups are not possible, the surgeon should be sent a control radiograph of the relevant joint. This will enable the surgeon to detect any complications at an early stage.

**6.4** CERAVER<sup>®</sup> acetabular components are compatible with MRI examinations. They may nonetheless create artefacts that should be taken into account for a perfect interpretation of these examinations. However, a period of 48 hours is recommended between implant placement and an MRI examination to avoid any movement of the device.

### 7- COMPLICATIONS

The complications of cemented dual mobility cups are as follows:

- Dislocation;
- Loosening;
- Osteolysis;
- Heterotopic ossification;
- Perioperative acetabular fracture;
- Traumatic fracture.

Complications related to the acetabular reconstruction cage:

- Fracture of the acetabular reconstruction cage or associated screws.

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