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BIPOLAR HEAD

Hemiarthroplasty Implant & surgical technique

AN EXTENDED RANGE (40 TO 56 MM) FOR ANY INDICATION OF HEMIARTHROPLASTY

EASY TO USE WITH :

- 22,2 mm diameter femoral head combined with stainless steel cups smaller than 48 mm
- 28 mm diameter femoral head combined with 48 mm stainless steel cups and greater

THE PERIPHERAL GROOVE HOLDS THE LOCKING RING AND PREVENTS ANY DISASSEMBLY WITHOUT USING PLIERS

CERAVER - Les Laboratoires Osteal Médical
69, rue de la Belle Etoile
95957 ROISSY CDG CEDEX - FRANCE
TEL. : +33 (0)1 48 63 88 63 - FAX : +33 (0)1 48 63 88 99
www.ceraver.com
contact@ceraver.com

BT_H29-02_FRGB_0721
GMED ISO 13485

SURGICAL TECHNIQUE

1 CEPHALIC CALIBRATION

Bipolar head gauge 40-46 mm
reference 24046

Bipolar head gauge 48-56 mm
reference 24056

Evaluate the bipolar trial cup size using the bipolar head gauges

2 TRIAL FEMORAL HEAD SIZE'S CHOICE*

Ø22,2 mm for sizes < 48 mm
Ø28 mm for sizes ≥ 48 mm

3 REDUCING WITH THE TRIAL BIPOLAR CUP

Correspondence table of trial bipolar cups

Ext diam. (mm)	Int diam. (mm)	Reference
ø 40	ø 22	2840
ø 42	ø 22	2842
ø 44	ø 22	2844
ø 46	ø 22	2846
ø 48	ø 28	2848
ø 50	ø 28	2850
ø 52	ø 28	2852
ø 54	ø 28	2854
ø 56	ø 28	2856

Trial assembly with the trial bipolar cup to pick the implant's final size

*Identical assembly with CERAFIT/OSTEAL trial instrumentation or CERAFIT/OSTEAL stem

4 FEMORAL HEAD ASSEMBLY

Manually clip the selected femoral head inside the {locking ring/liner} group

5 ASSEMBLY WITH THE STAINLESS STEEL CUP


Clip the previous assembly inside the stainless steel cup while checking the mobility of the femoral head and the locking ring

6 IMPACTION OF THE BIPOLAR HEAD

Impact the bipolar head on the selected femoral stem

IMPLANT

INSTRUMENTATION

CERAVER BIPOLAR HEAD

Reference	Size (mm)	Reference	Size (mm)
2808	40/22,2 mm	2012	48/28 mm
2809	42/22,2 mm	2013	50/28 mm
2810	44/22,2 mm	2014	52/28 mm
2811	46/22,2 mm	2015	54/28 mm
		2016	56/28 mm

TRAY REFERENCE 10436

OPTIONAL : FEMORAL HEAD REPLACEMENT

1' STAINLESS STEEL CUP DISASSEMBLY

Disassemble the {liner + locking ring + femoral head} group from the stainless steel cup using the bipolar pliers by closing the ring.

2' LINER DISASSEMBLY

Disassemble the locking ring out of the liner by opening the pliers and remove the femoral head.

3' FEMORAL HEAD REPLACEMENT

Replace the femoral head and start again from the surgical technique's step 4.