● NOTICE D'INSTRUCTIONS

Insert PE double mobilité

Français

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Dual mobility UHMWPE inserts

English

Instructions for use for the following implants:

- Dual mobility UHMWPE inserts (ref.10029; 10030; 10031; 10032; 10033; 10034; 10035; 10036; 10037)

Graphic symbols used for labelling these devices:

Symbol	Symbol title	Symbol	Symbol title
	Manufacturer		Do not re-use
	Expiration date		Device terminally sterilized by irradiation
	Batch code		Indicates a medical device having been submitted to a sterilisation process
	Catalogue reference of the device		N.B.: Not all the information required for use of the device is indicated on the label, it is necessary to consult the instruction manual

Prior to utilisation of CERAVER® hip joint prostheses, surgeons should familiarise themselves with the following recommendations:

1- GENERAL POINTS

1.1 The manufacturer and the supplier of the CERAVER® hip prosthesis provide a warranty for the design and material quality of the product. A prosthesis which has already been used must not be implanted again. Instructions for assembling the different components must be respected, as defined in the commercial documentation. It is strictly forbidden to use CERAVER® hip replacement components together with implants of another origin.

1.2 The following factors are of prime importance for the success of an operation using a CERAVER® hip prosthesis:

- Careful selection of a suitable patient for the operation in question.
- Exact surgical technique and proper use of special CERAVER® instruments.
- Extremely strict asepsis.

1.3 The surgeon is responsible for any complications that may result from an erroneous indication, incorrect surgical technique or inadequate aseptic precautions. Neither the manufacturer nor the supplier of the CERAVER® implants can be held liable for these complications.

1.4 The surgeon must be completely familiar with all the aspects and problems of hip replacement surgery. The surgical procedures for the implantation of a CERAVER® prosthetic hip can be best learned at a clinic with many years' expertise in the use of these implants.

1.5 It should be noted that for patients receiving a hip joint prosthesis that their weight and activity level may affect the longevity of the prosthesis.

2- DESCRIPTION OF IMPLANTS

CERAVER® cementless cup components are manufactured using the latest technology. Their mechanical and functional characteristics can only be guaranteed if the surgical technique and the present instructions are strictly observed.

- The inserts are made of ultra-high molecular weight polyethylene (UHMWPE) (ISO 5834-2 standards).

The angular movements afforded by CERAVER® double mobility cups are in compliance with NF EN ISO 21335.

3- INDICATIONS

- A joint prosthesis should be considered only after all other surgical methods of treatment have been carefully weighed up against each other and none have been judged to be more appropriate.

- Even the most successfully implanted artificial joint is inferior to a sound natural joint. On the other hand, an artificial joint can be a highly beneficial substitute for a deformed or diseased joint, because it eliminates pain and restores the restoration of good mobility and weight-bearing capacity.

- Every artificial joint is always subject to wear, which still remains a major problem that cannot be avoided. An initially stable prosthesis may become loose in the course of time. Wear and loosening of a prosthesis may make revision surgery necessary.

- Infection of a prosthesis will usually necessitate its removal or replacement in certain conditions and any such infection may leave debilitating residual effects.

Planned applications:

Hip arthroplasty consists of fully or partially replacing a hip joint suffering from various affections that cannot be treated through other therapies or methods. It aims to restore all or part of the patient's mobility and stability while minimising pain or even eliminating it.

The indications in the event of use of a double mobility PE insert in combination with the CERAVER® cementless Double Mobility HAC cup are as follows:

- Deterioration of the joint caused by primary osteoarthritis, inflammatory joint disease or developmental dysplasia of the hip;
- Traumatic fracture of the head or neck of the femur;
- Avascular necrosis of the femoral head;
- Risk of dislocation associated with muscle weakness or a dysplastic problem;
- Patients with a physical or mental deficiency that might increase their risk of hip dislocation;
- Unstable total hip replacement revision.

Based on the above statements, a joint prosthesis is indicated in cases where the following conditions are fulfilled:

- The selected patient's joint is anatomically and structurally suited for the reception and implantation of a prosthesis;
- Patients must have been informed by the surgeon of the risks associated with the implantation of a joint prosthesis and have declared their consent to the operation.

The following circumstances call for particular caution:

- Overweight patients;
- Excessive loading to be expected as a result of work and sport;
- Epilepsy or other factors favouring repeated accidents with increased risk of fracture;
- Severe osteoporosis and/or osteomalacia;
- Past history of still potentially active infectious diseases with potential joint manifestations;
- Severe deformation of the affected joint, which may render fixation of the implant more difficult;
- Mental inability of patients to understand and follow the surgeon's instructions;
- Drug or alcohol addition;
- Metabolic disorders or general medical treatments, leading to a progressive deterioration of bone support of the prosthesis;
- Other joint affections (knee, ankle, etc.);
- Patient's overall immune status weakened;
- Unbalanced diet.

4- CONTRAINDICATIONS

The following conditions are generally accepted as contraindications for the implantation of a joint prosthesis:

- Allergy to the implant materials of which the prosthetic components are made;
- Patients with a native cup with a diameter of less than 46 mm or greater than 62 mm;
- Active infection of the joint or surrounding area;
- Severe neuroarthropathy;
- Loss of bone structure proximally and distally of the joint, which precludes proper anchorage of the prosthesis;
- Significant muscle loss making joint prosthesis stability impossible, neuromuscular damage or vascular insufficiency of the limb affected, making the surgery unjustifiable;
- Charcot joint disease;
- Patients aged less than 50 years old and suffering from joint damage that could probably be treated successfully with a conservative treatment;
- Pregnancy;
- Tumor.

5- IMPORTANT WARNINGS

5.1 Preoperative warnings

- The operation should be planned carefully. On the basis of the radiological templates and pre-operative planning, the surgeon will choose the appropriate implant.
- The special CERAVER® templates and ancillary instrumentation should be used to prepare the fitting and setting of the artificial components. A description of the operative technique can never be exhaustive or contain all the risks and complications that should be considered.
- The cementless Double Mobility HAC cups must be used with a CERAVER® Double Mobility PE insert of the same diameter.

5.2 Intraoperative warnings

- Joint prostheses have to be handled with care. Even slight scratches or impact dents on the articular surfaces may contribute to excessive wear and initiate undesirable complications.
- To choose the right-sized implant: only CERAVER® templates, ancillary instrumentation and dummy implants must be used.
- Exposed articular surfaces must not come into contact with any metallic object (e.g. instruments, table tops, etc.).
- Assembly of the insert and corresponding head must be performed in accordance with surgical technique and using CERAVER® instruments.

6- POSTOPERATIVE PRECAUTIONS

6.1 All patients with an articular CERAVER® prosthesis must be postoperatively monitored by an orthopaedic surgeon. Follow-up examinations are essential to detect signs of wear or prosthesis loosening before any resulting detrimental effects. Early detection of an impending complication allows the surgeon to initiate timely and effective countermeasures. Above all, a revision surgery performed at the right time has a much better chance of success than a later surgical intervention, when advanced loosening has already resulted in serious resorption of the bone bed.

6.2 Patients must be instructed to inform their surgeon of the slightest changes in the joint operated on.

6.3 If annual clinical check-ups are not possible, the surgeon should be sent a control radiograph of the relevant joint. This will enable the surgeon to detect any complications at an early stage.

7- COMPLICATIONS

The complications in the event of use of a double mobility PE insert in combination with the CERAVER® cementless Double Mobility HAC cup are as follows:

Complications related to the cups:

- Loosening;
- Dislocation;
- Peri-prosthetic fracture;
- Peri-prosthetic osteolysis;
- Traumatic fracture;
- Wear.

Complications related to the surgical procedure:

- Superficial and deep infections;
- Haematoma;
- Claudication;
- Walking discomfort;
- Phlebitis;
- Deep vein thrombosis;
- Superficial necrosis of the wound;
- Groin pain.

8- COMPATIBLE MEDICAL DEVICES

Ancillaries:

- Dedicated instrumentation.
- CERAVER® stainless steel and alumina femoral heads with a diameter of 28 mm.
- CERAVER® cementless Double Mobility HAC cups.

9- STORAGE AND HANDLING OF PROTHESES

- Prosthetic components should be stored in their unopened original containers.
- Joint prostheses must not be subjected to treatment or modification.
- Prostheses that are damaged or have undergone a treatment that is inappropriate or not authorised by the manufacturer must not be implanted.
- Do not use the device after the expiration date. Contact your CERAVER® representative to return the product.

10- PACKAGING AND STERILISATION

The various components are packed individually in double packaging and are supplied sterile. All CERAVER® hip implants are sterilised by gamma radiation at a minimum dose of 25 kGy. The UHMWPE components are packaged under argon gas. It is important to check that the sterility protection is intact before use.

11- RESTERILISATION AND REUSE

Resterilisation and/or reuse of components of a CERAVER® hip prosthesis is strictly prohibited under penalty of exposing the concerned patient to contamination that can cause infection.

12- DISPOSAL AFTER USE

If a device is removed for revision surgery, please dispose of the components according to the centre's medical waste disposal procedure. This does not apply if the revision procedure is due to a Medical Device Vigilance case.

13- NOTE

The CERAVER® Double Mobility PE insert range obtained the EC mark in 2004.

Federal law in the USA and CANADA requires that prosthetic components be purchased by a pharmacist. For further information concerning the use of these components, please contact the manufacturer or its foreign subsidiaries and distributors.

CERAVER - Les Laboratoires Osteal Médical
69 rue de la Belle Etoile
CS 54263 - Gonesse - 95957 Roissy CDG Cedex - France
Tél. : 33 (0)1 48 63 88 63 - Fax: 33 (0)1 48 63 88 99
www.ceraver.com

Date of last revision of this document: March 2015

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