# RALLONGE TIBIALE

- **RALLONGE TIBIALE CERAVER**  
Notice d'instructions
- **CERAVER TIBIAL EXTENSION**  
Instructions for use
- **TIBIASCHAFT CERAVER**  
Anleitung
- **EXTENSIÓN DE LA TIBIA CERAVER**  
Manual de instrucciones
- **PROLUNGA TIBIALE CERAVER**  
Istruzioni per l'uso
- **ТИБИАЛЬНЫЙ УДЛИНЯЮЩИЙ СТЕРЖЕНЬ  
CERAVER**  
Инструкция по эксплуатации

LES LABORATOIRES OSTEAL MEDICAL

STERILE R

0459

## • RALLONGE TIBIALE CERAVER

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**• CERAVER TIBIAL EXTENSION**

Instructions for use

**Instructions for using the following component:**

Tibial extension	5080	5081	5082	5083	5084	5085	5086	5087	5088	5089
------------------	------	------	------	------	------	------	------	------	------	------

Symbol	Symbol title	Symbol	Symbol title
	Manufacturer		Do not reuse
	Expiry date		Please note that the special precautions for use should be checked for any important safety information such as warnings or precautions to be taken
	Batch code		
	Sterilized by irradiation		CE marking followed by the code for the notified body (0459, LNE/G-MED)
	Catalogue reference		

## • CERAVER TIBIAL EXTENSION - Instructions for use

Before using a CERAVER tibial extension, the attention of the surgeon is drawn to the following recommendations:

### 1 - GENERAL POINTS

- A. The tibial extension is an optional component of a CERAVER knee prosthesis and cannot be used on its own. Its use is left up to the assessment of the surgeon according to the characteristics of the patient (bone quality, age, physical activity, weight).  
The use of the CERAVER tibial extension with components of other origins is strictly prohibited.  
The manufacturer and supplier of a CERAVER tibial extension guarantee the quality of the manufacturing process and material used. A CERAVER tibial extension already used must not be implanted again.
- B. Several criteria are key to the success of a procedure to implant a CERAVER knee prosthesis:
  - Appropriate patient selection for the procedure;
  - Correct implant selection;
  - Observance of the operative technique and correct use of the surgical instruments specially designed by CERAVER;
  - Application of strict aseptic technique.
- C. The surgeon is responsible for any complications that may arise from erroneous indication and/or use of the instruments associated with the operative technique and/or of poor asepsis; neither the manufacturer nor the supplier of the CERAVER products can be blamed for these complications.
- D. The surgeon must be aware of all problems posed by implanting a knee prosthesis. The surgical procedure for implantation of a CERAVER total knee prosthesis can be best learned at a clinic or training center with surgeons who have acquired expertise in the technique.
- E. **The surgeon must inform patients receiving prosthetic knee joints of the fact that their weight and activity levels may influence the durability and performance of the prosthesis.**

### 2 - DESCRIPTION OF THE CERAVER EXTENSION

The tibial extension manufactured in titanium alloy (TiAl6V4 ELI - Standard ISO 5832-3) and is common to the HERMES PS, HERMES PS REVISION, HiFit PS and CERAGYR\* ranges.

The tibial extensions have a useful length of 50 mm. Their diameter increases by 1 mm steps from diameter 10 mm to diameter 19 mm.

Reference	5080	5081	5082	5083	5084	5085	5086	5087	5088	5089
Diameter (mm)	10	11	12	13	14	15	16	17	18	19

These tibial extensions can just as well be cemented or not cemented, but do not have any coating. The decision to cement extensions is left up to the surgeon's assessment.

### 3 - INDICATIONS

Refer to the instructions for use of the implant associated with the tibial extension.

The tibial extension is recommended in cases of severe deformations and insufficient bone quality to improve the anchoring of the tibial base.

#### 4 - EXPECTED PERFORMANCES

Refer to the instructions for use of the implant associated with the tibial extension.

#### 5 - CONTRAINDICATIONS

Refer to the instructions for use of the implant associated with the tibial extension.

#### 6 - WARNINGS

Refer to the instructions for use of the implant associated with the tibial extension.

The following factors can be extremely important for the success of the operation:

##### A. Instructions

Tibial extensions cannot be used apart from HERMES PS, HERMES PS REVISION, CERAGYR\* or HiFit PS tibial baseplates.

Extensions must be selected correctly to obtain good stability of the tibial baseplate. The degree of success in implanting the prosthesis is related to the correct selection of the diameter of the extension.

##### B. Handling the medical device

The exposed implant joint surfaces must not be placed in contact with any metallic objects or other hard objects that could potentially damage the surface.

##### C. Assembly and compatibility of HiFit PS and HERMES PS components

All the tibial extensions are compatible with all the HiFit PS and HERMES PS tibial baseplates.

The tibial extensions are optional components for HiFit PS and HERMES PS prostheses.

The keel of the HiFit PS and HERMES PS tibial baseplates has a threading plug (Ref. 6297) that must be removed before combining a tibial extension with the tibial baseplate.

Tibial extensions are assembled, in the case of HiFit PS and HERMES PS tibial baseplates, by a conical sleeve associated with an assembly screw (Ref. 3051):

- Assemble the extension in the baseplate without knocking it;
- Insert the assembly screw into the baseplate until it is screwed into the extension;
- Hit the extension at least 3 times with a 500 g hammer on a surface as rigid as possible while preventing metal-metal contact by interposing a compress folded 3 times;
- Finish the screwing operation.

##### CERAGYR\*

All the tibial extensions are compatible with all the CERAGYR\* tibial baseplates.

The tibial extensions are optional components for the CERAGYR\* prosthesis.

The keel of each CERAGYR\* tibial baseplate has a threading plug (Ref. 6299) that must be removed before combining an extension with the tibial baseplate.

Tibial extensions are assembled, in the case of CERAGYR\* tibial baseplates, by a conical sleeve: when the extension is in the tibial baseplate, hit the extension at least 3 times with a 500 g hammer on a surface as rigid as possible while preventing metal-metal contact by interposing a compress folded 3 times.

### **HERMES PS REVISION**

All the tibial extensions are compatible with all the HERMES PS REVISION tibial baseplates and all the offset adapters.

It is essential to associate an extension with HERMES PS REVISION prostheses (ref.: 5241, 5242, 5243, 5244, 5341, 5342, 5343, 5441, 5442, 5443, 5444):

- Assemble the extension in the baseplate without knocking it;
- Insert the assembly screw into the baseplate until it is screwed into the extension;
- Hit the extension at least 3 times with a 500 g hammer on a surface as rigid as possible while preventing metal-metal contact by interposing a compress folded 3 times;
- Finish the screwing operation.

An **offset adapter** (ref. 5806, 5807, 5808, 5809) can be inserted between the extension and the tibial baseplate. The tibial extension is then assembled with the offset adapter by conical sleeve:

- Once the offset adapter is positioned in relation to the tibial baseplate, insert the assembly screw;
- Position the offset adapter, as defined during the operation, in the baseplate without knocking it;
- Insert the assembly screw into the baseplate until it is screwed into the offset adapter;
- Place the extension in the offset adapter;
- Hit the extension at least 3 times with a 500 g hammer on a surface as rigid as possible while preventing metal-metal contact by interposing a compress folded 3 times;
- Finish the screwing operation.

## **7 - PRECAUTIONS**

Refer to the instructions for use of the implant associated with the tibial extension.

## **8 - ADVERSE EFFECTS**

Refer to the instructions for use of the implant associated with the tibial extension.

## **9 - PACKAGING AND STORAGE**

- A. The tibial extensions are stored in their original packaging.
- B. The various components are single-use, they are unit packaged in double bags and provided sterile.
- C. Do not use the device after the expiry date. Contact your CERAVER representative to return the product.

## **10 - STERILIZATION**

- A. All CERAVER knee implants are sterilized using gamma radiation at a minimum dose of 25 kGy; UHMWPE components are argon-conditioned to avoid polymer oxidation during sterilization.
- B. Before use, it is important to check the integrity of the sterility protection.

## **11 - RESTERILIZATION AND REUSE**

- A. Resterilization and/or reuse of any of the components of a CERAVER knee replacement is strictly forbidden because of the risk of contamination that could lead to infection.
- B. Resterilization and/or reuse of any of the components of a CERAVER knee replacement may compromise the device's mechanical performance and biological safety. **Never reuse surgical implants.** If an implant is removed after it has been implanted, it must be disposed of. Although it may appear intact, it may exhibit minor faults and signs of internal stresses that may undermine its integrity. Only new implants should be used and their configuration must not be modified before use.

## 12 - DISPOSAL AFTER USE

If a device is removed for follow-up replacement, dispose of the components according to the centre's medical waste disposal procedure. This does not apply if the follow-up procedure is due to a Medical Device Vigilance case.

## 13 - NOTE

The Ceraver knee prosthesis extension was granted CE Marking in 1998.

Federal law in the USA and CANADA requires that prosthetic components be purchased by a pharmacist. For further information concerning the use of these medical devices, please contact the manufacturer, its overseas subsidiaries or distributors.

\* CERAGYR: Knee prosthesis with mobile bearing tibial insert not approved in Canada

**CERAVER® - Les Laboratoires Osteal Médical**  
69 rue de la Belle Etoile  
CS. 54263 - 95957 Roissy CDG Cedex - France  
Tel.: 33 (0)1 48 63 88 63 - Fax: 33 (0)1 48 63 88 99  
[WWW.CERAVER.COM](http://WWW.CERAVER.COM)

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